

Solubility Data For Common International Colours

Colour	Distilled Water %	Propylene Glycol %	Ethanol %	Glycerine %
Tartrazine	31	06	Traces	18
Sunset Yellow	17	02	Traces	20
Ponceau 4R	30	10	Traces	0.50
Carmoisine	08	02.50	Traces	02.50
Amaranth	18	01.80	Traces	01.50
Brilliant Blue	28	12	15	20
Erythrosine	08	24	20	20
Indigo Carmine	01.50	0.70	Traces	1
Allura Red AC	22.0	0.50	Traces	3
Quinoline Yellow	14	14	Traces	-
Chocolate Brown HT	20	15	Traces	05
Black PN	08	01	Traces	0.50
Patent Blue V	06	02	Traces	03.50
Green S	15	01.50	Traces	01.50

This above data is obtained by dissolving the colour in hot solvent and allowing the resultant to stand for 24 hrs at room temperature